

MTP18-mitochondrial protein 18 kDa Gene

Reference: PTXBC001608

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	mitochondrial fission protein MTP18; MTP18; HSPC242; mitochondrial fission process protein 1; mitochondrial 18 kDa protein; mitochondrial protein 18 kDa; mitochondrial fission process 1		

For research use only.

GENE DESCRIPTION	mitochondrial protein 18 kDa
GENE ID	51537
ORIGIN SPECIES	Human
ACCESSIONS	BC001608
VECTOR	pDONR223
NCBI SYMBOL	MTP18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagagccgagccgctggggcgcagagcgcgatctaccgggacacgtgggtgcgatactggctatgccaatgaggtggcgagccttccgctctctgtgccagcggcggtgtgtgctgagctatggcgtggccagctcctacgtgctggcgatgccattgacaaag
gcaagaaggctggagaggtgccagccctgaagcagccgcagcggcagggtgacctgctgtgtggacacctttgatggcaggctctagcctctgtggccattccggcttcaccatcaaccgcgtgtgtgctgcctctctatgtcctggcactgccaccctggccctgg
ctgtccgcaagtggaccaccaccgcctgggctgtgacctcccatcattaccacccatgacagctcgtgtgatttctcctctggaactcagcctgcgcaagcttaccacagtgagggaagcccagctcctctga