



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

WDR34-WD repeat domain 34 Gene

Reference: PTXBC001614

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

DIC5; FAP133; SRTD11; bA216B9.3; WD repeat-containing protein 34; WD repeat domain 34

For research use only.

GENE DESCRIPTION	WD repeat domain 34
GENE ID	89891
ORIGIN SPECIES	Human
ACCESSIONS	BC001614
VECTOR	pDONR223
NCBI SYMBOL	WDR34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr