

RNASET2-ribonuclease T2 Gene

Reference: PTXBC001660

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNASE6PL; bA514O12.3; ribonuclease T2; ribonuclease 6		

For research use only.

GENE DESCRIPTION	ribonuclease T2
GENE ID	8635
ORIGIN SPECIES	Human
ACCESSIONS	BC001660
VECTOR	pDONR223
NCBI SYMBOL	RNASET2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccctgcagccctgcgcggggccctgctgggctgcctctgcctggcgttgccttgcctggcgtgcggaagcgcctgcgtgacaacctgagtggaataaactaatatggttcagcactggcctgagacagatgcgagaaaattcaaacgactgtagagacctccgat
tactggacaatacatggactatggccgataaaagtgaaggatgtaataagatcgtgcccctcaatttagaagagattaaggatctttgccagaaatgagggcactggtcctgacgttaattactcgtttcccaatcgagccgcttctggaagcatgagtggaagcatgggacctgcg
ccgcccagtggtgctgctcaactcccaagaagtaactttggcagaagcctggaactctacaggagctggacccaacagtgctctctaaattggggataaaacatccatcaattactacaggttgagattttaagatgccctggccagatgatggagtgatacccaaatcc
agtgcctccaccaagccagatgaggaagtacagacaattggcagatagaactgtgcctcaactaagcaagaccagcagctgcaaaactgcaccgagccgggggagcagcgtcccccaagcaggaagtctgctggcgaatggggcccgagagccgggtctgagagctctg
gaagatggccagctcttctatccccacctaataaagaccaagcattga