

**COX16-COX16 cytochrome c oxidase assembly homolog (S. cerevisiae) Gene**

Reference: PTXBC001702

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX16, cytochrome c oxidase assembly homolog; cytochrome c oxidase assembly protein COX16 homolog, mitochondrial; C14orf112; HSPC203; cytochrome c oxidase assembly factor		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	COX16 cytochrome c oxidase assembly homolog (S. cerevisiae)
<b>GENE ID</b>	51241
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC001702
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	COX16
<b>LENGTH</b>	Proteogenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggttcacccgcggtgatcgtcctttcgaagaacaagactctcggctatggagtcctcatgttggctgattgttggagttcctttgctctcgtgagtttctcaatccgatgatgctgtgaagagtaaaatggatcctgagcttgaagagataaaatcttttag  
agtcggaatatgagaaaatcaaagactccaagtgtgatgactggaagaatattcaggaccagccttgggaagatcctgacctccaaggaagaatccagaagccttaagactaagacaacttga