



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ANKRD10-ankyrin repeat domain 10 Gene

Reference: PTXBC001727

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat domain-containing protein 10; ankyrin repeat domain 10		
For research use only.			

GENE DESCRIPTION	ankyrin repeat domain 10
GENE ID	55608
ORIGIN SPECIES	Human
ACCESSIONS	BC001727
VECTOR	pDONR223
NCBI SYMBOL	ANKRD10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcgcgcgcgggagcggcgccgcggcgtagagcggcgcttcacgcagggagctcctctcgtcctgttcccgctgcaccgcgcctgcgcgcaggggacctgcgcacgctctgcctcgtcgcagcagacacccacgccacacgtgcctctcaggactcctctatgctcgtgcacgccctgctacggccgcgcacgttcggcagcttctagtcagttgtgtgagacgggagccacactcaacctctccaccacagctacgcgcagacgccagccacatgcagccttctggggacatcctcagtcctgtctgtcgtgattcaagcagagacccaattaa caaacccgattgtgagggtgaacctcccatcacaagcagctcgtcctcagccctgaatgcacatgcgcctctgtgcgggaatggcctcgaactgcagctgtattctacagagcagtgatcattctgtctccagacagtgatcattctgtatcctgaaccttaatttaattcagtgctcgaacccagcgttcgtctacagcagacatgtcctcaaacccagcgttctcccaagcgtgcgcctgttctcctgaagctcagcagctgtcatcctgaaccttctataa