

TFPT-TCF3 (E2A) fusion partner (in childhood Leukemia) Gene

Reference: PTXBC001728

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FB1; INO80F; amida; TCF3 fusion partner; INO80 complex subunit F; TCF3 (E2A) fusion partner (in childhood Leukemia); amida, partner of the E2A		

For research use only.

GENE DESCRIPTION	TCF3 (E2A) fusion partner (in childhood Leukemia)
GENE ID	29844
ORIGIN SPECIES	Human
ACCESSIONS	BC001728
VECTOR	pDONR223
NCBI SYMBOL	TFPT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagccgtggccttggaggattctcagcgcgcaggtcagagttggcgttgcctccctatttgggtgccacatcctggagagcagctggagacggaagtggagttgtgtcaggtggtctggcggtcagggctcgggagcagatgaagaggaagagcagccgg
ggctggcgcgccgagcgggaattaaatgcagaagtaccagcactaggtcggcgtcggggagatcgaacgagcgggtcctgaacaggtccatcaggtgagaggataactcggaggtcagcaggaacgaggttctcatgagagtgctgactcct
acggggatgactaccggccagccagttcaccattgtgtgaggtgagggcagccagggcacggaatgccccaccacaggaatgcggagaatgagcctccagagaagagacactgtccccccgaaggactcctgcacccccagaacccggcagccagccccgg
gaggggcccagtgggcggaagagcggcagtgccacgggatggagcgcgagcaggaatgcgtgactccagagctggccccggtgcagattaaggttgaggagactttgcttgaagcagatgagggcctgattccagttgggtttctggggtccagacaactgtgcc
ctacccgacctggccagccagcctctgactga