

TOMM7-translocase of outer mitochondrial membrane 7 homolog (yeast) Gene

Reference: PTXBC001732

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TOM7; mitochondrial import receptor subunit TOM7 homolog; 6.2 kd protein; translocase of outer membrane 7 kDa subunit homolog; translocase of outer mitochondrial membrane 7 homolog; translocase of outer mitochondrial membrane 7		

For research use only.

GENE DESCRIPTION	translocase of outer mitochondrial membrane 7 homolog (yeast)
GENE ID	54543
ORIGIN SPECIES	Human
ACCESSIONS	BC001732
VECTOR	pDONR223
NCBI SYMBOL	TOMM7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaagctgagcaagaggccaagcagagactacagcagctcttcaaggggagccagtttgcattcgctggggctttatccctcttgatgtacctgggatttaagaggggtgcagatcccgaatgcctgaaccaactgtttgagcctactttggggataa