



Reference: PTXBC001746

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BPTP-4; HEPTP; LC-PTP; LPTP; PTPNI; tyrosine-protein phosphatase non-receptor type 7; dual specificity phosphatase 1; hematopoietic protein-tyrosine phosphatase; protein-tyrosine phosphatase LC-PTP; protein-tyrosine phosphatase, nonreceptor-type, stress i		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, non-receptor type 7
GENE ID	5778
ORIGIN SPECIES	Human
ACCESSIONS	BC001746
VECTOR	pDONR223
NCBI SYMBOL	PTPN7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]