



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C15orf39-chromosome 15 open reading frame 39 Gene

Reference: PTXBC001762

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

uncharacterized protein C15orf39; chromosome 15 open reading frame 39

For research use only.

GENE DESCRIPTION	chromosome 15 open reading frame 39
GENE ID	56905
ORIGIN SPECIES	Human
ACCESSIONS	BC001762
VECTOR	pDONR223
NCBI SYMBOL	C15orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagcagcccccattggccacaggagctgtgtgaagcaagtgcagtgacttggacctgaggctggagctgtgtctcttcagagccccactgtggccagagatgagccagagagcgcttagcccttggtcagaagtaccggcccccgaagctcaggaagccaggcagaagaagccac
caaccctctggccgggagaagaagcagagcagctgctggagagcagctctgtgtgtgtctcccttaccctgtaacacagctgccagcggccactgggccccactgaagcggccctctccgcatctgtctggagaccgcctgtgtctaafggctctgctgtgccactgtgggccaagctct