

RPL34-ribosomal protein L34 Gene

Reference: PTXBC001773

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L34; 60S ribosomal protein L34; leukemia-associated protein; ribosomal protein L34		
For research use only.			

GENE DESCRIPTION	ribosomal protein L34
GENE ID	6164
ORIGIN SPECIES	Human
ACCESSIONS	BC001773
VECTOR	pDONR223
NCBI SYMBOL	RPL34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtccagcgtttgacataccgacgttaggctttcctacaatacagcctcctacaaaactaggctgtccgaaccctggtaatagaattgtttacctttataccaagaaggttgggaaagcaccaaaatctgcatgtggtgtgtgccaggcagacttcgaggggttcgtgctgtaagacctaa
agttcttatgagattgtccaaaacagaacatgtcagcaggcctatgtgttccatgtgtgctaaatgtgttcgtgacaggatcaagcgtgtttcctatcaggagcagaaaaatcgttgaaggtgtgaaggcacaagcacagagtcagaaagctaaataa