



Reference: PTXBC001794

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAG3A; CTG3a; NBL; NUMB-R; NUMBLIKE; NUMBR; TNRC23; numb-like protein; numb homolog-like; numb-related protein; NUMB like, endocytic adaptor protein		

For research use only.

GENE DESCRIPTION	numb homolog (Drosophila)-like
GENE ID	9253
ORIGIN SPECIES	Human
ACCESSIONS	BC001794
VECTOR	pDONR223
NCBI SYMBOL	NUMBL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr