

GSG1-germ cell associated 1 Gene

Reference: PTXBC001796

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	germ cell-specific gene 1 protein; testicular secretory protein Li 20; germ cell associated 1		

For research use only.

GENE DESCRIPTION	germ cell associated 1
GENE ID	83445
ORIGIN SPECIES	Human
ACCESSIONS	BC001796
VECTOR	pDONR223
NCBI SYMBOL	GSG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctctgaaggccttctctggccagcggacactcctatctgccatctcagcatgctatcagcttctccacaacatccctgctcagcaactactggtttgtgggcacacagaaggtgcccaagccctgtgcgagaaaggtctggcagccaagtgccttgacatgccagtgtcc
ctggatggagataccaacacatcccccaggaggtggtacaatacaactgggagactggggatgaccggttctcctccggagctccggagtgcatgtgctatctgtgaggaaactgtggagaaccaggggagaggtgccgaagtctcatgaacttacaccaccagccaagag
aggtgagaaaggactactggaattgccacgttgcaaggcccatgtcaccctcctccgatttgagggaagcggttgatggagaaggcttccctccctccctccctggggccttgggcaaaaatcctatggttatccctgggaacgcagatcacctacatcggaactcaattcatcag
cttctcctgctactaacaagactgtactactcctgggaaccctgcctgtgggctcaaaactgagcgcttctgctgttctcctgtcctgtcaggtctcctggggatgggtggccacatgatgtattcacaagcttccaagcgactgtcaactgggtccagaagactggagaccacatgtttgg
aattatggctgggcttctacatggcctggctccttcacctgctgcatggcgtcggctgtcaccacctcaacacgttacaccagatgggtgctggagtcaagtgcaagcatag