

SNX6-sorting nexin 6 Gene

Reference: PTXBC001798

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MSTP010; TFAF2; sorting nexin-6; TRAF4-associated factor 2; tumor necrosis factor receptor-associated factor 4 (TRAF4)-associated factor 2; sorting nexin 6		

For research use only.

GENE DESCRIPTION	sorting nexin 6
GENE ID	58533
ORIGIN SPECIES	Human
ACCESSIONS	BC001798
VECTOR	pDONR223
NCBI SYMBOL	SNX6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgaaggagaattcacaagatgaacagggaactggaagctgaatattggcaatattcaagaagacagttgcgatgcatgaagtgtctgtgtcgtgtggcagcacatctatttgagaagagattaaattccatgtcttctggaataatacaagatttgagtgctgcgaggaataataaaaaagagaacttgaaagactctttaaatacatggttaaatcagcagatggagtaaatcgtttcaggagtaaggatgtagatgatttcttgagcagcagcaacatttcttggaaatcataaccgagtttaaggatgcatctgctaaatctgatagaatgacaagatcccaaaaagtgcgcagatgattacaatagaattggttcttattatgcttttaggaactcaggattctacagatataatgcaagtttttctcaagtttcagaactgttcgataaaacaagaaaaatagaagcagcagtgctgtctgatgaagacctcaaaccttctgactctttaaataattacttaagagaatctcaagctgctaaggatctctgtatcgaaggtctaggtcactagtggaattatgaaatgctaataaagcactggataaagcagagcaaaaaataaagatgttctacaggccgaacttcccaacaattatgtgtcagaaattgaaaaaatatctgagctgcaaaacaagaacttatagattttaagacaagaagattgctgcattcagaaaaatttagtggaaactggcagagttagaactgaagcatgcaaaaggtaattctacagttgctgcagaactgcctggcagtggttaattggagacacataa