

CDKN2D-cyclin-dependent kinase inhibitor 2D (p19, inhibits CDK4) Gene

Reference: PTXBC001822

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INK4D; p19; p19-INK4D; cyclin-dependent kinase 4 inhibitor D; CDK inhibitor p19INK4d; cell cycle inhibitor, Nur77 associating protein; cyclin-dependent kinase 4 inhibitor D p19; cyclin-dependent kinase inhibitor 2D (p19, inhibits CDK4); inhibitor of cycli		

For research use only.

GENE DESCRIPTION	cyclin-dependent kinase inhibitor 2D (p19, inhibits CDK4)
GENE ID	1032
ORIGIN SPECIES	Human
ACCESSIONS	BC001822
VECTOR	pDONR223
NCBI SYMBOL	CDKN2D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtgaggaggttcgcgcggcgaccggctgagtgggcgggcgccggcgagctgcaggaggtgcggcctctgcaccgcgagctgggtcatccgcgcctcaaccgcttcggcaagacggcgtgcaggctcatgattggcagcaccgccatcgccctgg
agctgtgaagcaagggtccagcccaatgtccaggacacctccggtaccagtcagtcacatgacgcagcccgactggattcctggacacccfagaggtcctagtgagacacgggctgatgtcaacgtgctgatggcaccggggcacttcaatccatctggcagttcaagaggt
cacactgtgtgtcagctttctggcagctgaatctgatctccatcgaggagacccaggggtctcacaccttggagctggcactgcagaggggctcaggacctctgtggacatcctgcagggccacatggtgccccgtgtga