

TRIM48-tripartite motif-containing 48 Gene

Reference: PTXBC001862

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNF101; tripartite motif-containing protein 48; RING finger protein 101; tripartite motif containing 48		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 48
GENE ID	79097
ORIGIN SPECIES	Human
ACCESSIONS	BC001862
VECTOR	pDONR223
NCBI SYMBOL	TRIM48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattctggaatctgcaagctctccagagggaactcacctgccccatctgcatgaactactcatagaccgggtcaccatagactgtgggcacagctttgagggccctgtttctacctaactggcaagacatcccaattcttactcagtgctttgaatgcataaagacaatacagcagaga
aacctcaaaactaacattcgattgaagaagatggcttcccttgcagaaaagccagctctgctattctgagctctgaggagcaaatgtgtggcattcacaggagacaaagaagatgttctgtgaagtgacaggagcctgctgtgtgtgtctccagctctcaggagcaccgggtatc
acagacactgtcccgtgagtggtgctgctgaggaacactgggagaagctgttaagaaaatgcagctttatgggaaaaagcttgtgaaaatcagagaaacctgaatgtggaaccaccagaatcagccactggaaggcttttgagacatattatcaggagtgagtcctgtgctgca
catgccccagccttgaaatctagcgtcagggcagggccatcactggactgaggacaggctcaaccaattctga