

## COPS3-COP9 constitutive photomorphogenic homolog subunit 3 (Arabidopsis) Gene

Reference: PTXBC001891

|                |                                                                                                                                                                                                                        |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CSN3; SGN3; COP9 signalosome complex subunit 3; COP9 complex subunit 3; COP9 constitutive photomorphogenic homolog subunit 3; JAB1-containing signalosome subunit 3; signalosome subunit 3; COP9 signalosome subunit 3 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | COP9 constitutive photomorphogenic homolog subunit 3 (Arabidopsis)                               |
| <b>GENE ID</b>          | 8533                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC001891                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | COPS3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcgtctgccctggagcagttcgtgaacagtgctccgacagctctcagctcaaggcgaatgacacagctttgtgaactgatcaacaagagtgagggaactcttgcgaagaacttaccatctggacactgtgctcggggctctggatgtacaagaacactcctgggcgtctgtgttt  
tgtttgtgaagttttctatgccagtggtctgacttcgaaacgctattctcaggttgcagctcttcacagcactgtgaatggggagcacttcgatatgcaacagacacttttgcgggctttgccatcagtaacaaatgcacttgggaagaaaacagccctgcgaggaattggcctccta  
agcaagccatagacaagatgcagatgaatacaaacagctgacctcaatacatgctgctcctgcccagctttgtttgctagcaaaatgctttaagcctgccctccatcttgacgtggatgatggaatctgtaagagaatggagcctgatgcaaaacacttttattgttactattattgg  
agggatgatttactggcgtgaagaactttgaagagctctctacttttgaacaggctataactcctgcccagctcagtcagtcataatcgttggaatcatataaaagtatatatttagtgcctttgatattacttggcaagtacaacagctaccataatacatctcaaatgtgggtagatt  
cattaagcccttagcaatgcataccagagtagcacaagtgtattcaaccaacaacccctcagaactccgaaacctggtgaataagcacagtgaaaccttcactcgcgataacaacatggggctggggaagcaatgcttgcactctttataagaagaatttcagaggctacaagaagacc  
tttttaactctatcattacaagatggcaagctgctgagcttgcagctcaggaggcagagaatactgtctgcacatgatagaagatgggtgagattttgcaagtattaaccagaaggacggtatggtcagtttccatgataacctgaaaaatataataaccagccatgctcataacat  
tgatcaggagatgctgaagtgcttgagctggatgagcggctgaagccatggaccaggagatcacagtgaacctcagttgtacaaaagatgaggctcacagaagatgattcaggaaacaacacatcagttattcttga