

IL10RB-interleukin 10 receptor, beta Gene

Reference: PTXBC001903

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDW210B; CRF2-4; CRFB4; D21S58; D21S66; IL-10R2; interleukin-10 receptor subunit beta; IL-10 receptor subunit beta; IL-10R subunit 2; IL-10R subunit beta; IL-10RB; cytokine receptor class-II CRF2-4; cytokine receptor class-II member 4; cytokine receptor f		

For research use only.

GENE DESCRIPTION	interleukin 10 receptor, beta
GENE ID	3588
ORIGIN SPECIES	Human
ACCESSIONS	BC001903
VECTOR	pDONR223
NCBI SYMBOL	IL10RB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtggagcctgggagctggcgtgggtgctgctggtgtcagcattgggaatgtaccaccccgaaatgtcagaatgaattctgttaattcaagaacattctacagtgaggagtcacctgctttgccaaagggaacctgactttcacagctcagtaacgttatagatattc
caagataaatgcatgaatactaccttgacggaatgtgatttctcaagctcttccaaagtatgggtgaccacaccttgagagtcagggtgaatttgcagatgagcaatcagactgggttaacatcacctctgtctgttgatgacacattatggacccctggaaatgcaagtagaagtactgtg
attctttacataatgctgttcttagccctaaatgagaatgaatacgaacttggactatgaagaatgtgtataactcatggactataatgtgcaatactggaaaacggtactgatgaaaagtccaattactccccagatgactttgaggtctcagaaaacctggagccatggacaactatg
tgttcaagttcaggggtttctctctgacggacaagaactggggaatggatgagcctgtctgtgagcaacaacccatgacgaacggtccctctggtggtgcccgtatcctcatggcctcgtctcatggtctgctgctggcactcctcggctctcgtctgctgtgtgtgctgttaca
agaagacaaagtacgcttctccctaggaattctctccacagcactgaaagagtttttggccatcctcatcataacacactctgttttctcttccattgtcggatgagaatgatgttttgacaagcgaagtgtcattgcagaagactctgagagcggcaagcagaatcctgtgtgacagct
gcagcctcgggaccccgctgggcagggtggcccaagctag