

MTX1-metaxin 1 Gene

Reference: PTXBC001906

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MTX; MTXN; metaxin-1; mitochondrial outer membrane import complex protein 1; metaxin 1		

For research use only.

GENE DESCRIPTION	metaxin 1
GENE ID	4580
ORIGIN SPECIES	Human
ACCESSIONS	BC001906
VECTOR	pDONR223
NCBI SYMBOL	MTX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgcccatggagctgttctgctgctcagggggctgggggctgcccgtcagtgacctggacagcctggccgtgctgacattgctgctccactgaaggtacacaaagatcagcaaccctggcagagcccttcaggaactctgctgcccctcgaccagtcacat
ggagaggtcatctcagttccacacaaagatcagccaccccttcgaaaagaggtacatacttttggatagacaccaaactacgtggaagtgaacccggaagtggatgacagagctatgcccttccctcaacttctctgctggccgcatgcagcggcagatggaacggctaca
gctgctgactgggggacacagggcctgaggacgaggaagagctggagaaggagctgtaccgagagggctcgggagtgctgacccctgctctcagcgcctgggctcctcaaaagtcttcttggagatgccctgctccttggacgcttctgcttgcagctacttggccctgctgctgag
gcaaaagctgccagtgaggagctgcaggtccactgctggtgggtgcacaacctctgtgcttaltgtaccacattctcagttctctacttcccctgggatggagctgaggtaccacgcacgacagacagcagggccagagactgaggaggagccataccggcgccggaaccag
atcctatctgtgctgagcaggtgagcagccatgggtggctacgctgtcagcggcattgtctccatccagcgggcaacgcctgctcgggccccagggcaccctggcctgagtgaggatgaagaggatga