

C16orf14-chromosome 16 open reading frame 14 Gene

Reference: PTXBC001912

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf14; FAM195A; c349E10.1; MAPK regulated corepressor interacting protein 2; MAPK regulated co-repressor interacting protein 2; family with sequence similarity 195, member A; protein FAM195A		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 14
GENE ID	84331
ORIGIN SPECIES	Human
ACCESSIONS	BC001912
VECTOR	pDONR223
NCBI SYMBOL	C16orf14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtacaccatcaccaaggggcccagcaagctggtcgcgagcgccgcacaggtcccacgcagcagcaggtggaggccggctggcgagctcctgaaatgccgcagcccgccgcgacctgcagccccgcggcgagcccttgcgcagccgcccggaccctggc
ccctgtcgagtcagggccaaggcctgtgttcaatcgtgtgaatggccggcgcccccctccacgtcccaccccttcgaggggacccaggagacctacacagtgccccacgaggagaatgtccgctttgtgtccgaagcctggcagcaggtgcaacagcagctggatgggtggccagc
cgggtgaggcgccgccaaggccctgtgcagtagctggagaggaccccaatccccggctgcagaactttgtgccattgacctagacagagtggtggcgagcagttcctggcgagaatcaccagctgttcctag