



Reference: PTXBC001915

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C8orf35; INT10; integrator complex subunit 10		
For research use only.			

GENE DESCRIPTION	integrator complex subunit 10
GENE ID	55174
ORIGIN SPECIES	Human
ACCESSIONS	BC001915
VECTOR	pDONR223
NCBI SYMBOL	INTS10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]