

C1orf116-chromosome 1 open reading frame 116 Gene

Reference: PTXBC001943

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	specifically androgen-regulated gene protein; specifically androgen-regulated protein; chromosome 1 open reading frame 116		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 116
GENE ID	79098
ORIGIN SPECIES	Human
ACCESSIONS	BC001943
VECTOR	pDONR223
NCBI SYMBOL	C1orf116
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagatgtgagggcagatgaaggaatgttctggaagggtgcctttaaggctgctgagaatttctggactgtgtcctgatggacgcagaccatcaagcccagaatttctga