

PNRC2-proline-rich nuclear receptor coactivator 2 Gene

Reference: PTXBC001959

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	proline-rich nuclear receptor coactivator 2; proline rich nuclear receptor coactivator 2		
For research use only.			

GENE DESCRIPTION	proline-rich nuclear receptor coactivator 2
GENE ID	55629
ORIGIN SPECIES	Human
ACCESSIONS	BC001959
VECTOR	pDONR223
NCBI SYMBOL	PNRC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtggtggagaggtataacattccagccctcaatctagaaatgttagtaagaaccaacaacagcttaacagacagaagaccaaggaacagaattccagatgaagattgttcataagaaaaaagaagaggacatggtataactcatcagcagctgcctggcaggccatgcaa
aatgggggaagaacaaaaatttccaataatcaaaagttggaattctagcttatcaggtccaggttacttttaaatctcaagctaatcagaactatgctggtgccaaatttagtgagccgccatcaccaagtgcttctccaaccacaaagccactgggtccctgtttcccttaactcttcagat
aaggaaataatgacatttcaacttaaaccttaactaaagtacaggtataa