



Reference: PTXBC001961

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SNAT7; amino acid transporter; solute carrier family 38 member 7		
For research use only.			

GENE DESCRIPTION	solute carrier family 38, member 7
GENE ID	55238
ORIGIN SPECIES	Human
ACCESSIONS	BC001961
VECTOR	pDONR223
NCBI SYMBOL	SLC38A7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]