

TAF1D-TATA box binding protein (TBP)-associated factor, RNA polymerase I, D, 41kDa Gene

Reference: PTXBC001972

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JOSD3; RAFI41; TAF(I)41; TAFI41; TATA box-binding protein-associated factor RNA polymerase I subunit D; Josephin domain containing 3; RNA polymerase I-specific TBP-associated factor 41 kDa; TATA box binding protein (TBP)-associated factor, RNA polymerase		

For research use only.

GENE DESCRIPTION	TATA box binding protein (TBP)-associated factor, RNA polymerase I, D, 41kDa
GENE ID	79101
ORIGIN SPECIES	Human
ACCESSIONS	BC001972
VECTOR	pDONR223
NCBI SYMBOL	TAF1D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggataaatcaggaatagattctcttgaccatgtgacatctgatctgtggaacttgcaaatcgaagtataactctctgatagcagcttatttaaaactcagtgatccctactcacctaaaggggagaaaagaaccccatcgaaaattgttctgacacctgaaagtgttcacgcaagtga
ttcatcaagtgaactcattttgaaccaataccattgactataaaagctattttgaaagattcaagaacaggaagagagatataaaaaaagaaaagaggaggtaccagccaacaggaagaccacggggaagaccagaaggaggaatcctatactactactaataagataagaag
aaacaatttagaagcagaggatctggcttccatttttagaatcagagaatgaaaaaacgcaccttggagaaaaatttaacgtttgagcaagctgtgtcgaagaggatttttaactatattgaaaaactgaagtatgaacaccacctgaaagaatcattgaagcaaatgaatgttggaagattt
agaaaatgaagatttgacagtcgtagatacaaattttgatgatgatggatccattctctattgaggagtcacagcagaggatgaggatgcaacacatcttgaagataacgaatgtgatcaaatggcaggggatgttcatagttaagttctgaattccctgtaagactgagtgatact
agaagaagaggatattactgaagaagctgcttctctaaaaagagagctacaaaagccaaaatactggacagagagggcctgaaaatgtga