

C6orf106-chromosome 6 open reading frame 106 Gene

Reference: PTXBC002328

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf106; FP852; dJ391O22.4; chromosome 6 open reading frame 106		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 106
GENE ID	64771
ORIGIN SPECIES	Human
ACCESSIONS	BC002328
VECTOR	pDONR223
NCBI SYMBOL	C6orf106
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcatggactgtagactggaccggagctgatgcagaagttcagctccctggcaccacgacaaggagctgctcatctccgagttccagaggtgctcggcttcagctcaatcctgcccgttgcccttctctggacatgaccaactggaactacaagcagcaattggc
gcctattatgactttgagagcccaaacatcagtgccctctatgctctttggaagatgaccatagagagaaggaggagtcataactccggatactcagttgttaaaaacatggcgatccagaattctgggagagcctggcctccagggtttgtcttaaatatgctggggagacc
aattggacatgtaacatggtgatggtgagatcgtagagcccaagagatgacagatgtagcgtccagatgtagcagcccagcagagcaggaatgatacaggagacagtgccgagtgctactgctacaggactctactatggagatgcatctgggtgatctcagttggagggtggg
tggacttttaggagtaacgagcagctgctcctttgaaacggagtcaacacacagccgcatgtaaggtagaaggaaactcaaccctttgctctcccaaaagaaccgacaatcagatgaaaacaactaaaagaccctggggctccgagttcgactcagcagcaaaacacat
gggctcctgctctgacacatgggctcctgctcctgaccaaaactgagcaagaccagaatagactgacagaactctgtaaatctgctcccagcagtcacgaaaacaactatcagtagtacttacagtaagggtccatgggcttaccctcggccagcttaa