

C11orf79-chromosome 11 open reading frame 79 Gene

Reference: PTXBC002331

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf79; PGL2; SDH5; succinate dehydrogenase assembly factor 2, mitochondrial; SDH assembly factor 2; hSDH5; succinate dehydrogenase subunit 5, mitochondrial; succinate dehydrogenase complex assembly factor 2		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 79
GENE ID	54949
ORIGIN SPECIES	Human
ACCESSIONS	BC002331
VECTOR	pDONR223
NCBI SYMBOL	C11orf79
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtgtctacagtggtctcgaactcgtcgtgatgcttgcctctgtcaaggcacagcctattgtctctttgtcagtgtagacatcattcagacgcttctacagaggtgacagcccaacagattccaaaaggacatgattgaaatcccttgcctccatggcaggagagaactgatgaatcca
tagaaacaaaagagcccgctgctctatgagagcagaagagggaatgttgaaaactgcattctcttagtcttttgcataaagaacatctgcagacatgacagaaaagcagctgaacctctatgaccgctgattaacgagcctagtagatgagggatattactactggccacaga
agctaaaccagccccagaaatatttgaaatgaagtcattgcccctgctgagagactttgctaaaaacaaaacaaagacagagactgcgtgccccagatcttgagtaccttttgaaaagccacgttga