

PSMD9-proteasome (prosome, macropain) 26S subunit, non-ATPase, 9 Gene

Reference: PTXBC002383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Rpn4; p27; 26S proteasome non-ATPase regulatory subunit 9; 26S proteasome regulatory subunit p27; homolog of rat Bridge 1; proteasome (prosome, macropain) 26S subunit, non-ATPase, 9; proteasome 26S non-ATPase regulatory subunit 9; proteasome 26S subunit,		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) 26S subunit, non-ATPase, 9
GENE ID	5715
ORIGIN SPECIES	Human
ACCESSIONS	BC002383
VECTOR	pDONR223
NCBI SYMBOL	PSMD9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgacgaggaagcgaggcagagcggaggtcctcgcaggccggcgtcgtgactgtcagcgacgtccaggagctgatcggcgcaaggaggagatagaagcgcatcaaggccaactatgacgtgctgaaagccaaaaggcattgggatgaacgagccgctgtgga
ctgtgagggtctacccccgttcagacgtggacgtgtaccaagtccgacccgagggcacaacatcatatgcctgcagaatgatcacaaggcagtgatgaagcaggtggaggaggccctgcaccagctgcacgctcgcgacaaggagaagcagcccggtgacatggctgagggccac
aaagagcccatgagccgcaactgggtcagagtgagagccaggccctccacgggcttcgcaaaagtgaacagcatcagccccggctcccgagccagcgcgggtctgcaagtggatgatgagattgtggagttcggctctgtgaacaccaagaactccagtcactgcataaca
ttggcagtggtgacgacagtgagggggaagccctgaatgtgacagtgatccgaggggggaaaaacaccagcttagactgttccaacacgctgggcaggaaaaggactgctgggctgcaacattattcctctgcaaatgatga