

## NAP1L1-nucleosome assembly protein 1-like 1 Gene

Reference: PTXBC002387

|                |                                                                                                                                                             |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | NAP1; NAP1L; NRP; nucleosome assembly protein 1-like 1; HSP22-like protein interacting protein; NAP-1-related protein; nucleosome assembly protein 1 like 1 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | nucleosome assembly protein 1-like 1                                                              |
| <b>GENE ID</b>          | 4673                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC002387                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | NAP1L1                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcagacattgacaacaaagaacagctgaacttgatcaagatttggatgatgttgaagaagtagaagaagaggaaactggggaagaacaaactcaaacgacgtcagctaacgttcagatgatgcaaaatcctcagattctgcagccctcaagaagacttgatggtctgtagaa  
acaccaacaggatgacattgaaagcctgctaggtgagttaaagacgagtgatgctctcaaaacctgcaagttaaatgtgcacagatagaagccaaattctatgaggaggttcagcttgaaaggaggtatgctgttctctatcagcctctattgataagcgatttgaattattaatgcaat  
ttatgaacctacggaagaagaatgtgaatggaaccagatgaagaagatgagatttggaggaattgaaagaaagcccaagattgaagatgagaaagagatgaagaaaaagaccccaaggaattcctgaatttgggttaactgttttaagaatgttgacttgctcagtgatggtt  
caggaaacacgatgaacctattctgaagcacttgaagatataaagtgaagtctcagatgctggccagcctatgagtttctctagaatttcaacttgaacccaatgaatattttacaatgaagtgtgacaagacatacaggatgaggtcagaaccagatgattctgatcccttttttgatg  
gaccagaattatgggtgtacaggtgcccagatagattggaagaaaggaagaatgtcatttgaactatataagaagaagcagaacacaaagggacgtgggacagttcgtactgtgactaaacagtttccaatgactctttcttaactttttgcccctcctgaagttcctgagatggag  
atctggatgatgctgaagctatccttgcagacttcgaattgtcactttttacgtgagcgtataatcccaagatcagtggttatatttactggagaagctattgaagatgatgatgattatgatgaagaaggtgaagaagcgatgaggaaggggaagaagaaggatgaggaa  
aatgatccagactatgacccaagaagatcaaaacccagcagagtgcaagcagcagtg