

HMOX2-heme oxygenase (decycling) 2 Gene

Reference: PTXBC002396

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	heme oxygenase 2; heme oxygenase (decycling) 2		
For research use only.			

GENE DESCRIPTION	heme oxygenase (decycling) 2
GENE ID	3163
ORIGIN SPECIES	Human
ACCESSIONS	BC002396
VECTOR	pDONR223
NCBI SYMBOL	HMOX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagcgggaagtggaaacctcagaggggtagacgagtcagaaaaaagaactctggggccctagaaaaaggagaacaaatgagaatggctgacctctcgagctctgaaggaaggacaaaggaagcacacgacgggcagaaaaacaccagttgtcaaggactcttgaa
aggcaacattaagaaggagctgtttaagctggccaccagggcactttacttcacatactcagccctcaggaggaatggagcgcaacaaggacatccagcctttgcccctttgacttcccatggagctgcaccggaaggaggcgtgaccaaggacatggagatattcttggtgaaa
actgggagagcaggtgcagtgcccaaggctgcccagaagtagtgagcgggacccatcacatagggcagaaacgagcgggagctactgtggcccatgcatacacccgctacatgggggatctctcgggggcccaggtgctgaagaagggtggccagcagcactgaaactccc
cagcacaggggaaggacccagttctactgtttgagaatgtggacaatggccagcagttcaagcagctctacgggcccagatgaacgcccgtgacctaacaatgaagaccaagagagagatcgtggagggagccaacaaggcttttgagtataacatgcagatattcaatgaactgg
accaggccggctccacactggccagagagacctggaggtgggtccctgtacacgatgggaaggagacatgcgtaaatgccctttctacgtgctgaacaagcaaaaggtgccctggaggcagcagctgtcccttcgaacagatggtgtgctgaggaagcccagcctcca
gttcacctggccgctggtgtggccctagctgctgagactctggcctgttactacatgtga