



Reference: PTXBC002408

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASC-1; TARBP-B; 60S ribosomal protein L3; HIV-1 TAR RNA-binding protein B; ribosomal protein L3		

For research use only.

GENE DESCRIPTION	ribosomal protein L3
GENE ID	6122
ORIGIN SPECIES	Human
ACCESSIONS	BC002408
VECTOR	pDONR223
NCBI SYMBOL	RPL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]