

HTATIP2-HIV-1 Tat interactive protein 2, 30kDa Gene

Reference: PTXBC002439

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	oxidoreductase HTATIP2; CC3; SDR44U1; 30 kDa HIV-1 TAT-interacting protein; HIV-1 TAT-interactive protein 2; HIV-1 Tat interactive protein 2, 30kDa; Tat-interacting protein (30kD); short chain dehydrogenase/reductase family 44U, member 1; HIV-1 Tat intera		

For research use only.

GENE DESCRIPTION	HIV-1 Tat interactive protein 2, 30kDa
GENE ID	10553
ORIGIN SPECIES	Human
ACCESSIONS	BC002439
VECTOR	pDONR223
NCBI SYMBOL	HTATIP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaacagaagccctgtcgaagcttcgggaagactcaggatgcagaataatccgtcttattttggcgccagcgggagaaccggcagagtgtcttaagggaatcctggagcagggcctgtttccaaagtcacgtcattggccggaggaagctcaccttcgacgaggaagcttataaaatgtgaatcaagaatgtgtgacttggaaatgtgactacgcctctgccttcaaggtcatgatgttgattctgttgccctgggtaccaccagagggaagctggggcggagggaattgttcgtgttgaccgagattatgtgctgaagtctgcagagctggcaaaagctggagggtgcaaacattcaacttgcctccttaaggagctgataaatcaagcaatttttatatctacaagttaagggaagtagaagccaaggttgaagaattaaattgatcgttactctgtatttagccctggagttctgttatgtataggcaagaatctgccagggtgaatggctgttagaaattcttggctctaccagactctggccagggggcattctgtgcctgtgtgaccgtgttagagcaatgctgaacaatgtgttgagaccaagagacaagcagatgaactgtggagaacaaggccatccatgacctggggaagcgcagctctctcaagccatga