

DNAJB6-DnaJ (Hsp40) homolog, subfamily B, member 6 Gene

Reference: PTXBC002446

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DJ4; DnaJ; HHDJ1; HSJ-2; HSJ2; LGMD1D; LGMD1E; MRJ; MSJ-1; dnaJ homolog subfamily B member 6; DnaJ (Hsp40) homolog, subfamily B, member 6; DnaJ-like 2 protein; heat shock protein J2; DnaJ heat shock protein family (Hsp40) member B6		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily B, member 6
GENE ID	10049
ORIGIN SPECIES	Human
ACCESSIONS	BC002446
VECTOR	pDONR223
NCBI SYMBOL	DNAJB6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggattactatgaagttctagggcgtgcagagacatgcctcaccgagggatataaaaggcatatcggaaactggcactgaagtggcatccagataaaatcctgagaataaagaagaagcagagagaaaattcaagcaagtagcggagggcatatgaagtgcctggatgtaag
aaacgggacatctatgacaaatattggcaagaaggatataatgggtggaggaggagtggaagtcattttgacagtcatttgaatttggcttcacatccgtaaccagatgatgtctcagggaatttttgggtggaagggaaccattttcatttgactctttgaagacccttttgaggactctttg
ggaatcgaagggtcccgaggagcagaagccgaggggacggggcgtttttctcgttcagtggaatttccgtcttttgaagtggaatttctctttgatacaggatttactcatttgggtcactaggtcacgggggctcacttcattctctccacgtcatttgggttagtgcatgggc
aactcaaatcgatatacaactcaactaaatgggttaatggcagaaaaaacactacaagaagaattgtcgagaacgggtcaagaagaagtagaaggtgaagaagatggccagttaaagtccttaacaataaatgggtggccgacgacgatgccctcgtgaggagcgcacgagagggcc
agaacgccctgccagcccagcctgccggcctccgcccgaagccgcccggcctgcctcgtcgtgagacacgcgcctcactgtctctctgaggaggaggcgagcaggaccgacctcgggcacccgggcccctgggacccctcgcgtccgcagcagatgaagaaggtg
gcaagaggagaagcagaagcagagagaggagtcgaagaagaagaagtcgaccaaaggcaatcactag