

DYNLRB1-dynein, light chain, roadblock-type 1 Gene

Reference: PTXBC002481

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BITH; BLP; DNCL2A; DNLC2A; ROBLD1; dynein light chain roadblock-type 1; ROBL/LC7-like 1; bithoraxoid-like protein; dynein, cytoplasmic, light polypeptide 2A; dynein-associated protein Km23; roadblock domain-containing protein 1		

For research use only.

GENE DESCRIPTION	dynein, light chain, roadblock-type 1
GENE ID	83658
ORIGIN SPECIES	Human
ACCESSIONS	BC002481
VECTOR	pDONR223
NCBI SYMBOL	DYNLRB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaggtggaggagacactgaagcgactgcagagccagaaggagtgcaagggaatcctcgtcgaacacagaaggcattccatcaagagcaccatggacaacccccaccaccaccagatgccagcctcatgcacagcttcatcctgaaggcagggagcaccctgcgtgacatcgacccccagaacgatcctcctcgaattcgtccaagaaaatgaaattatggtgcaccagataaagactatttctgattgtgattcagaatccaaccaataa