

RPP25-ribonuclease P/MRP 25kDa subunit Gene

Reference: PTXBC002497

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ribonuclease P protein subunit p25; RNase P protein subunit p25; ribonuclease P 25kDa subunit; ribonuclease P/MRP 25kDa subunit; ribonuclease P/MRP subunit p25		

For research use only.

GENE DESCRIPTION	ribonuclease P/MRP 25kDa subunit
GENE ID	54913
ORIGIN SPECIES	Human
ACCESSIONS	BC002497
VECTOR	pDONR223
NCBI SYMBOL	RPP25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaactccgtaaggtgcgtccgaagaggcgccagcggggtgcggggccgaggaggcgccggggctccggcccttcgcagacctggcgccggcgcggtgcacatcggggtcaaggaggcagcaagatccggaacctgatggccttcgccaccgccagcatgg
cgagccagccacgcgcgcatcgtcttcagcggtgcggcgccgaccacaaaaccgtcacgtgcggcgagatcctcaagcgccgctggcgggcctgcaccaggtcacgcggctgcgtaccggagcgtaacgcgaggtgtggcagagcctccgctggggccacgcag
ggtcagacgcctggcgagccggcgtatgtcagcgtacttaagaacgtgcccggcctgcacatcctacttccaaggacgcgctggatccgcgacagcccggctaccagcccccgaatccccatcctgtccctcgtcccgccagccgcgacgctcaagaggagcctaggg
gaacccgcagctggagaaggctccggaagcgatcgcaacccgagccaggggttcggacgaggatcagacggcctga