

PDCD10-programmed cell death 10 Gene

Reference: PTXBC002506

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCM3; TFAR15; programmed cell death protein 10; TF-1 cell apoptosis-related protein 15; apoptosis-related protein 15; cerebral cavernous malformations 3 protein; programmed cell death 10		

For research use only.

GENE DESCRIPTION	programmed cell death 10
GENE ID	11235
ORIGIN SPECIES	Human
ACCESSIONS	BC002506
VECTOR	pDONR223
NCBI SYMBOL	PDCD10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggatgacaatggaagagatgaagaatgaagctgagaccacatcatgtttctatgccctctatgcagtcagtcctgtgttaatgagctagaacgagtaaatctgtctgcagccagacactgagagccgtttcatcaaggctgaaaaagaaaatccaggtctcacacaagac
atcattatgaaaatttagagaaaaaagcgtggaagttaacttcacggagtccttctctgtaggcagctgatgtagaagagtatatgattgaacgaccagagccagaattccaagacctaaacgaaaaggcagcagcacttaacaaattctcagtaagatccagatgagatcaatg
acagagtgaggttctgcagacaatcaaggatagctagtgcataaaagaactcttgatacagtgtaataatgtcttcaagaaatatcaatacagaacgcagggcacttgaacacaaaagaagaattgttaagtactccaaagtctcagtgatactctgaaaacgtattttaagatg
gcaaggcaataaatgtgttcgaagtgccaaccgactaattcatcaaaccaacttaacttcagacctcaaaactgtggcctga