

PIR-pirin (iron-binding nuclear protein) Gene

Reference: PTXBC002517

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pirin (iron-binding nuclear protein)		

For research use only.

GENE DESCRIPTION	pirin (iron-binding nuclear protein)
GENE ID	8544
ORIGIN SPECIES	Human
ACCESSIONS	BC002517
VECTOR	pDONR223
NCBI SYMBOL	PIR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggtcctccaagaagtactctcagtcagccgggagcagtcggaaggggttgagcaggggtccggagaagcattggcagaccgagttaaaaatctggatccgttttactgtttgatgaatttaaaggagtagaccaggaggatttctgatccacatcaggtttt
gaaacagtatctactctcctggaagggggcagcatggcccatgaagactctgtggacacactggtaaaatgaaccagggagatttgcagtggaatgactcggggccgggcatctgcacgtgagatgcttgcagaggagccagccatggcctacaactgtgggttaattgagg
agctcagagaagatgggtgagcctcagtcaccaggaactgaaaagtgaagaatccctaaaccagtaaggatggtgtgacagttgctgcatcttctggagaagccctgggaataaagtcgaaggtttacactcgcacaccaacttatatttggacttcaaatggacccaggagccaaaca
ttcccaactatccctaaaggggtggacaagcttcatctacagatatctggagatgtgtatattggcccgatgatgcacaacaaaaatagaacctcatcacagcagtcgttggagaaggtgacaggtgtccaggtggagaacaaggatcccaagagaagccacttcttaattgctgg
ggagccattaagagaaccagttatccaacatgtccatttgtgatgaacaccaatgaagagatttctcaagctattcttgatttcagaacgcaaaaaatgggttgaaggggccaaaacctggaaatcaagattgggaactag