

C20orf149-chromosome 20 open reading frame 149 Gene

Reference: PTXBC002531

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf149; dJ697K14.9; exdpf; pancreatic progenitor cell differentiation and proliferation factor; exocrine differentiation and proliferation factor; pancreatic progenitor cell differentiation and proliferation factor homolog		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 149
GENE ID	79144
ORIGIN SPECIES	Human
ACCESSIONS	BC002531
VECTOR	pDONR223
NCBI SYMBOL	C20orf149
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccatccccccagcggctcgtcgtggccacccagactactaccggcgccgctgggtccacttcagcaacagctcctgcagcagtagtgcggcgaggagccattccccccccaggctcccaaggctgacccgggtcattgggtggccagcttcttttc
gggaagtccaccctccggttcattggccacgggtgttgagtcgcagagcactcggaaactccccaggcctccagcagcatgaccgctgtggcctggctcgggacgccccgagggaagcagcccgccggtcagtcacagcagccagcgtgggccccgtcctga