

CHCHD7-coiled-coil-helix-coiled-coil-helix domain containing 7 Gene

Reference: PTXBC002546

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX23; coiled-coil-helix-coiled-coil-helix domain-containing protein 7; COX23 cytochrome c oxidase assembly homolog; coiled-coil-helix-coiled-coil-helix domain containing 7		

For research use only.

GENE DESCRIPTION	coiled-coil-helix-coiled-coil-helix domain containing 7
GENE ID	79145
ORIGIN SPECIES	Human
ACCESSIONS	BC002546
VECTOR	pDONR223
NCBI SYMBOL	CHCHD7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctcggtaacacagaggtgagagatcctgacataaatcctgtttgtcgggaatctgatgctccaccagatgtctggatgaaataactatgacagggaagggtgttcacttactcttgaggtacaaaaactgccggagattctggaattctatcgtgatgcagagaagaaagaacgg
agtgaagccatttatgcctacggcagcagaagagatgaaatcttgagagcagtgaggaaatgcctattga