



Reference: PTXBC002562

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARC41; p40-ARC; p41-ARC; actin-related protein 2/3 complex subunit 1B; ARP2/3 protein complex subunit p41; actin related protein 2/3 complex, subunit 1B, 41kDa; arp2/3 complex 41 kDa subunit; actin related protein 2/3 complex subunit 1B		

For research use only.

GENE DESCRIPTION	actin related protein 2/3 complex, subunit 1B, 41kDa
GENE ID	10095
ORIGIN SPECIES	Human
ACCESSIONS	BC002562
VECTOR	pDONR223
NCBI SYMBOL	ARPC1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr