

STAG3L2-stromal antigen 3-like 2 Gene

Reference: PTXBC002581

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGC131759; STAG3L1; STAG3L3; OTTHUMP00000209501; OTTHUMP00000209503; OTTHUMP00000209504; STAG3-like protein 2; stromal antigen 3-like protein 2; stromal antigen 3-like 2		

For research use only.

GENE DESCRIPTION	stromal antigen 3-like 2
GENE ID	442582
ORIGIN SPECIES	Human
ACCESSIONS	BC002581
VECTOR	pDONR223
NCBI SYMBOL	STAG3L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcaccgagaagtcgctgaagtcgctgaaggctctgaaagggctgtacgtaaccgggacgtgacgcacgcctggagctcttcactggccgcttcaaggactggatggtttccatgatcgtggacagagagtagactgtggcagtgaggccgtcagattactgatactatc
cttaagaacatggaaggggtgctgatggacgtggactgtgagagcgtctacccattgtgtag