

TMEM85-transmembrane protein 85 Gene

Reference: PTXBC002583

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM85; PIG17; ER membrane protein complex subunit 4; cell proliferation-inducing gene 17 protein; transmembrane protein 85		

For research use only.

GENE DESCRIPTION	transmembrane protein 85
GENE ID	51234
ORIGIN SPECIES	Human
ACCESSIONS	BC002583
VECTOR	pDONR223
NCBI SYMBOL	TMEM85
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggcccaagggggcctgggtgctaaccgaggcggcgttcaagtggccattgagctaagcgggcctggaggaggcagcaggggtcgaagtaccggggcagtgccaggagagactcgtctaccagtcgggtacttgacaagcaagtgcctgataccagcgtgcaa
gagacagaccggatcctgggtggagaagcgtcgtgggacatcgcttgggtccctcaaacagattcccatgaatctcttcacatgtacatggcaggcaatactatctccatcttccctactatgatgtgtgtatgatggcctggcgaccattcaggcacttatggccattcagccaagaa
tggagttcagtggtggaggactgctttgtgaacatgagaaagcagcgcctggcctatgtattgggtcttattacatcctctttaa