



Reference: PTXBC002641

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CK15; K15; K1CO; keratin, type I cytoskeletal 15; CK-15; cytokeratin 15; keratin 15, type I; keratin-15, basic; keratin-15, beta; type I cytoskeletal 15; keratin 15		

For research use only.

GENE DESCRIPTION	keratin 15
GENE ID	3866
ORIGIN SPECIES	Human
ACCESSIONS	BC002641
VECTOR	pDONR223
NCBI SYMBOL	KRT15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1