



Reference: PTXBC002666

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanylate-binding protein 1; GBP-1; GTP-binding protein 1; guanine nucleotide-binding protein 1; guanylate binding protein 1, interferon-inducible; guanylate binding protein 1, interferon-inducible, 67kDa; huGBP-1; interferon-induced guanylate-binding pro		

GENE DESCRIPTION	guanylate binding protein 1, interferon-inducible, 67kDa
GENE ID	2633
ORIGIN SPECIES	Human
ACCESSIONS	BC002666
VECTOR	pDONR223
NCBI SYMBOL	GBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1