

MED18-mediator complex subunit 18 Gene

Reference: PTXBC002694

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	p28b; mediator of RNA polymerase II transcription subunit 18; TRAP/mediator complex subunit p28b; mediator of RNA polymerase II transcription, subunit 18 homolog; mediator complex subunit 18		

For research use only.

GENE DESCRIPTION	mediator complex subunit 18
GENE ID	54797
ORIGIN SPECIES	Human
ACCESSIONS	BC002694
VECTOR	pDONR223
NCBI SYMBOL	MED18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcaccctccagtcacatgatgctgtcactggggcaccattaaatgatggagctgttgcagggaagtgttttagatcacagtttggaagcctccaccgccttcgtgtgtgtgacacatggaaacctgagacttcttgaccatgagatggfattcctccttaaggc
cagcaagccagccattgttctcaggcccgacgtctatggacagggcagggcaccctggcatctgcctacctgggacagccagaatgggagacaagaaccgcatgcctggfgcgaaactgcgtggacattgccacatctgagaacctcaccgactcttgatggaatgg
gctccgcatgacatgagttgttgtaaggacattgttccgtaaggcatcatgaagattatggtgtacaagattttccgcatcctgggtccagggaacacagacagcactgaggcctgtcactctcctatctcgtggaattaaagtgtgtagcaccgctgggacagatggtctct
gatgacatgaagaacttcgagaacagctaaaacctctggttcacctagagaaatagacccaagaggctcatgtga