

## AAK1-AP2 associated kinase 1 Gene

Reference: PTXBC002695

|                               |                                                                                       |                     |                |
|-------------------------------|---------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | AP2-associated protein kinase 1; adaptor-associated kinase 1; AP2 associated kinase 1 |                     |                |
| <b>For research use only.</b> |                                                                                       |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | AP2 associated kinase 1                                                                           |
| <b>GENE ID</b>          | 22848                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC002695                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | AAK1                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaagaagtttttgcactccggcgagagcagggcggtctgctccgctccagcggaggagggggcagcacctcgggcctgggcagtggtctacatcggaagagctctcggcatcggcgacagcaggtcacagtgagcaggtgttgccggaaggtggattgctatt  
gtatttctgtgaggaacaagcaatgggatgaaatgtgcttgaacgcatgtttgtcaacaatgagcatgatctccaggtgtgcaagagagaaatccagataatgaggatcttccaggcacaagaatattgtgggttacattgattctagtatcaacaacgtgagtagcgggtgatgtagga  
agtgtctattctgtaggactttttagaggtggccaggtgtgtaaacctgatgaaccagcgcctgcaaacaggtttacagagaatgaagtgctcagatattttgtgatacctgtgaagctgttgcctcctgcatcagtgcaaaactcctattatccaccgggacgtgaaggtgaaacatcc  
tcttgcagaccgagggcactatgtctgtgtgactttggaagcggccacaacaaatccagaatccacaactgagggagtgcaatgcagtagaagatgagattaagaatacacaaacgtgtcctatcgagcaccagaatgggtcaacctgtacagtggtgcaaaatcatcactacgaagga  
gacatttggcctctgtgattgtgtgtataaattgctacttcttgcatttgggaaagtcaggtggcaatttgtgtggaacttcacaattctgataattctgatatctcaagacatgcactgcctaataggtatattgttgaaccagacctgacaaaaaggccgatattaccaggt  
gtcctacttctcatttaagctactcaagaaagagtgcccaattccaaatgtacagaactctccattctgcaaaagcttctgaaccagtgaaagccagtgaggcagctgcacaaaagaccagccaaaggccagactgacagatccctccaccacagagactcaattgcacccggc  
agaggcctaaagctggcgagactcagcgaaccaggaaatcttccatccagccagcgtgacacccgggaagagggccactgttgcagccccacctcagcgtgcaggtatccagcaatcagcctggccttttagccaggtgttcccaacaaaacccaagccccccagccagc  
ctctgcccaaacctcagcccaagcagccacaggtctctccactccacagcagacgcttctactcaggccagggtgtgcccgtcagggccagccacacccagcaccagcagcagatacaataaaacttagtatgaaacttt