



Reference: PTXBC002701

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASCIZ; ZNF822; ATM interactor; ATM Interacting protein; ATM/ATR-Substrate Chk2-Interacting Zn ⁺⁺ -finger protein; ATM/ATR-substrate CHEK2-interacting zinc finger protein; ATM/ATR-substrate CHK2-interacting zinc finger protein; zinc finger protein 822		

[illegible]