

SNX5-sorting nexin 5 Gene

Reference: PTXBC002724

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sorting nexin-5; sorting nexin 5		

For research use only.

GENE DESCRIPTION	sorting nexin 5
GENE ID	27131
ORIGIN SPECIES	Human
ACCESSIONS	BC002724
VECTOR	pDONR223
NCBI SYMBOL	SNX5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcgggtcccgagttgctgcagcagcaggaggaggaccgcagcaagggtgaggtcctccagcctcagactcctggcgccggcgccctgcgagctccgggtcgtgcattcttcccttgcgcgtccataggccgtggctgtcgcgcgcacatccccagcgcgggaagcacc
gtgagacccggcgccctctccctggtcttctcactgaaggctgccccggggaatccctgtggtatgctccggattctattaggacagaaccagagacgcgggacctagccccagctcagcgccgggtgccgagcggacttggagagatgatctctggagaccccggaatgtctttta
aaacttctctcggcttcgtgtga