

DAP-death-associated protein Gene

Reference: PTXBC002726

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DAP-1; death-associated protein 1; death associated protein		
For research use only.			

GENE DESCRIPTION	death-associated protein
GENE ID	1611
ORIGIN SPECIES	Human
ACCESSIONS	BC002726
VECTOR	pDONR223
NCBI SYMBOL	DAP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcgcctccgaagggaactagagactaaagctggacacccgccgcgtgaaagctgggtggaatgcgaattgtgcagaaacacccacatacaggagacaccaaagaagagaaagacaaggatgaccaggaatgggaagccccagtcacctaaccactgtgttcac
tctggggatcagcccggtgacaaagattccccccggtgctgcaggtggtcaccagaagccgcatgcctccatggacaagcatcttcccaagaacccagcacatccagcagccacgcaagtga