



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

PLLP-plasma membrane proteolipid (plasmolipin) Gene

Reference: PTXBC002760

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PMLP; TM4SF11; plasma membrane proteolipid (plasmolipin); transmembrane 4 superfamily member 11 (plasmolipin)		

For research use only.

GENE DESCRIPTION	plasma membrane proteolipid (plasmolipin)
GENE ID	51090
ORIGIN SPECIES	Human
ACCESSIONS	BC002760
VECTOR	pDONR223
NCBI SYMBOL	PLLP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]