

UBE2G1-ubiquitin-conjugating enzyme E2G 1 (UBC7 homolog, yeast) Gene

Reference: PTXBC002775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E217K; UBC7; UBE2G; ubiquitin-conjugating enzyme E2 G1; E2 ubiquitin-conjugating enzyme G1; ubiquitin carrier protein G1; ubiquitin conjugating enzyme E2G 1; ubiquitin-conjugating enzyme E2G 1 (UBC7 homolog, C. elegans); ubiquitin-conjugating enzyme E2G 1		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2G 1 (UBC7 homolog, yeast)
GENE ID	7326
ORIGIN SPECIES	Human
ACCESSIONS	BC002775
VECTOR	pDONR223
NCBI SYMBOL	UBE2G1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggagctgcagtcggcactgctactgcgaagacagctggcagaactcaacaaatccagtggaaagctttctgcaggtttaatagatgacaatgatctctaccgatgggaagtccttattggccctccagatacattatgaagtggtgttttaaggctcatcttacttcccaa
aagattatccctccgacctcctaaatgaaattcattacagaaatctggcaccctaaatgttgataaaatggtgatgtgtgcatttcttctcatgagcctggggaagataagtatggttatgaaaagccagaggaaacgctggctccctatccacactgtggaacatcatgattagtgcatt
tctatgctggcagaccctaatggagactcacctgctaattgttgatgctgcgaagaatggagggaagatagaaatggagaatttaaagaaaagttgcccgcgtgtgtaagaaaaagccaagagactgcttttgagtga